



XL FDP	Latex Agglutination Slide Test
Application	For the detection and semiquantitation of cross linked fibrin degradation products in plasma
Sensitivity	≥ 200 ng/ml
Specificity	IgG class of Murine Monoclonal antibodies, specific to cross linked fibrin degradation products
Latex reagent	Uniform suspension of polystyrene latex particles coated with mouse monoclonal, Anti-D dimer antibody
PBS Buffer	For performing semiquantitation
Positive Control	Reactive with the reagent
Negative Control	Non reactive with the reagent
Test duration	Three minutes assay
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
XL FDP kit contains reagents and accessories	15 Tests	10650015
	60 Tests	10650060

