



Sensitive Thromboplastin Reagent with 3.2% Buffered Citrate Blood collection anticoagulant

Application	Combipack of Uniplastin [®] (PT) reagent and 3.2% buffered citrate with zwitterionic buffer for achieving precise, consistent results in PT determinations.
Reagent	<ul style="list-style-type: none"> • Ready to use Liquid stable, Calcified Thromboplastin reagent • 3.2% Buffered Tri sodium citrate
Source	Rabbit Brain
Normal Range	11 - 15 seconds
ISI value	Optimised at ~ 1.0 (1.0 -1.15)
Storage / Stability	2-8°C, 18 months

Presentation	Pack	Cat.No.
UNIPLASTIN PRECISION COMBIPACK with Uniplastin and 3.2% buffered citrate	5 ml (25 Tests)	10621025

