



	Chemiluminence assay for the detection of antibodies to Hepatitis C virus in human serum or plasma
Application	For detection of antibodies to Hepatitis C Virus in human serum or plasma.
Principle	Chemiluminence (CLIA)
Reagents	<ul style="list-style-type: none"> • Chemiluminence compatible microwells coated with and recombinant antigens for HCV • Controls for assay validation • Conjugate for antigen & antibody detection • Chemiluminescent substrate • Wash buffer
Storage / Stability	2-8°C, 12 months

Presentation	Pack	Cat.No.
Electra HCV Ab CLIA kit contains: <ul style="list-style-type: none"> • Microwell strips (3 x 8 wells/ pouch) • Reagents • Accessories 	96 Tests	407020096
	192 Tests	407020192

