



	<p>Fourth generation Chemiluminensence assay for the detection of antigens & antibodies to HIV 1&2 virus in human serum or plasma</p>
<p>Application</p>	<p>For detection of antibodies to HIV V1/2 & O subtypes and p24 antigen in human serum or plasma.</p>
<p>Principle</p>	<p>Chemiluminensence (CLIA)</p>
<p>Reagents</p>	<ul style="list-style-type: none"> • Chemiluminensence compatible microwells coated with agglutinating sera and recombinant antigens for both HIV-1 and HIV-2 • Controls for assay validation • Conjugate for antigen & antibody detection • Antibody Reagent: Biotinylated agglutinating sera for p24 antigen. • Chemiluminescent substrate • Wash buffer
<p>Storage / Stability</p>	<p>2-8°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>Electra HIV Ag/Ab 4.0 CLIA kit contains:</p> <ul style="list-style-type: none"> • Microwell strips (3 x 8 wells/ pouch) • Reagents • Accessories 	<p>96 Tests</p>	<p>407010096</p>
	<p>192 Tests</p>	<p>407010192</p>

