



	<p><b>Chemiluminence assay for the detection of antibodies to <i>Treponema Pallidum</i> in human serum or plasma</b></p>
<p><b>Application</b></p>	<p>For detection of total antibodies (i.e IgG, IgA &amp; IgM etc) to <i>Treponema Pallidum</i> in human serum or plasma.</p>
<p><b>Principle</b></p>	<p>Chemiluminence (CLIA)</p>
<p><b>Reagents</b></p>	<ul style="list-style-type: none"> <li>• Chemiluminence compatible microwells coated with recombinant 47Kd &amp; 17Kd antigens</li> <li>• Controls for assay validation</li> <li>• Conjugate for antibody detection</li> <li>• Chemiluminescent substrate</li> <li>• Wash buffer</li> </ul>
<p><b>Storage / Stability</b></p>	<p>2-8°C, 12 months</p>

Presentation	Pack	Cat.No.
<p><b>Electra Tp. Ab. CLIA</b> kit contains:</p> <ul style="list-style-type: none"> <li>• Microwell strips</li> </ul>	96 Tests	407040096
	192 Tests	407040192

