



<p style="text-align: center;">FIAcheck™ Anti-CCP</p>	<p>Fluorescence Immunoassay for quantitative determination of IgG class auto antibodies against cyclic citrullinated peptides in human serum and plasma</p>
<p>Application</p>	<p>FIAcheck Anti-CCP Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of IgG class auto antibodies against cyclic citrullinated peptides in human serum and plasma.</p>
<p>Principle</p>	<p>Rapid, lateral flow, quantitative immunofluorescence assay</p>
<p>Measuring range</p>	<p>8-500 U/mL</p>
<p>Reagents</p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> •Immobilized CCP antigen on test line •Goat anti chicken IgY in control line •Mixture of goat anti-human IgG antibody and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
<p>Test duration</p>	<p>15 mins</p>
<p>Storage / Stability</p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>FIAcheck Anti-CCP kit contains an individually pouched device, sample diluent, calibration QR code card and dilution vials</p>	<p>10 Tests</p>	<p>543010010</p>