



FIAcheck™ **CRP**

Fluorescence Immunoassay for quantitative determination of CRP in human serum, plasma or whole blood

Application	FIAcheck CRP Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of CRP in human serum, plasma or whole blood.
Principle	Rapid, lateral flow , quantitative immuno-fluorescence assay
Measuring range	0.5-200µg/mL
Reagents	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> •Immobilized CRP mouse monoclonal antibody1 on test line •Sheep anti chicken IgY in control line •Mixture of CRP mouse monoclonal antibody2 and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
Test duration	3 mins
Storage / Stability	4-30°C, 12 months

Presentation	Pack	Cat.No.
FIAcheck CRP kit contains an individually pouched device, sample diluent, calibration QR code card and dilution vials	25 Tests	537010025