



<p style="text-align: center;">FIAcheck™ LH</p>	<p>Fluorescence Immunoassay for quantitative determination of Luteinizing hormone (LH) in human serum</p>
<p>Application</p>	<p>FIAcheck LH Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of Luteinizing hormone (LH) in human serum.</p>
<p>Principle</p>	<p>Rapid, lateral flow, quantitative immunofluorescence assay</p>
<p>Measuring range</p>	<p>1.02-200mIU/mL</p>
<p>Reagents</p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> •Immobilized LH mouse monoclonal antibody1 on test line •Goat anti chicken IgY in control line •Mixture of LH mouse monoclonal antibody2 and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
<p>Test duration</p>	<p>15 mins</p>
<p>Storage / Stability</p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>FIAcheck LH kit contains an individually pouched device and calibration QR code card</p>	<p>25 Tests</p>	<p>535020025</p>