



<p style="text-align: center;">FIACheck™ PRL</p>	<p>Fluorescence Immunoassay for quantitative determination of Prolactin (PRL) in human serum</p>
<p>Application</p>	<p>FIACheck PRL Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of Prolactin (PRL) in human serum.</p>
<p>Principle</p>	<p>Rapid, lateral flow, quantitative immuno-fluorescence assay</p>
<p>Measuring range</p>	<p>25-4000µIU/mL</p>
<p>Reagents</p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> •Immobilized PRL mouse monoclonal antibody1 on test line •Sheep anti chicken IgY in control line •Mixture of PRL mouse monoclonal antibody2 and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
<p>Test duration</p>	<p>15 mins</p>
<p>Storage / Stability</p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>FIACheck PRL kit contains an individually pouched device and calibration QR code card</p>	<p>25 Tests</p>	<p>535030025</p>