



<p style="text-align: center;">FIACheck™ TSH</p>	<p>Fluorescence Immunoassay for quantitative determination of TSH (Thyroid Stimulating Hormone) in human serum and plasma</p>
<p>Application</p>	<p>FIACheck TSH Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of TSH (Thyroid Stimulating Hormone) in human serum and plasma.</p>
<p>Principle</p>	<p>Rapid, lateral flow, quantitative immuno-fluorescence assay</p>
<p>Measuring range</p>	<p>0.1-60µIU/mL</p>
<p>Reagents</p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> • Immobilized TSH monoclonal antibody 1 on testline • Sheep anti chicken IgY in control line • Mixture of monoclonal TSH mouse monoclonal antibody 2 and chicken IgY labeled with fluorescent microspheres on the binding pad <p>Sample diluent</p>
<p>Test duration</p>	<p>15 mins</p>
<p>Storage / Stability</p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p>FIACheck TSH kit contains an individually pouched device, sample diluent, calibration QR code card and dilution vials</p>	<p>25 Tests</p>	<p>53401025</p>