



<p><b>FIAcheck™</b> <b>TT3</b></p>	<p><b>Fluorescence Immunoassay for quantitative determination of TT3 (Total Triiodothyronine) in human serum and plasma.</b></p>
<p><b>Application</b></p>	<p>FIAcheck TT3 Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of TT3 (Total Triiodothyronine) in human serum and plasma.</p>
<p><b>Principle</b></p>	<p>Rapid, lateral flow, quantitative immuno-fluorescence assay</p>
<p><b>Measuring range</b></p>	<p>0.5-10nmol/L</p>
<p><b>Reagents</b></p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> <li>• Immobilized T3-BSA full antigen complex on test line</li> <li>• Sheep anti chicken IgY in control line</li> <li>• Mixture of T3 sheep monoclonal antibody and chicken IgY labeled with fluorescent microspheres on the binding pad</li> </ul> <p>Sample diluent</p>
<p><b>Test duration</b></p>	<p>15 mins</p>
<p><b>Storage / Stability</b></p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p><b>FIAcheck TT3</b> kit contains an individually pouched device and special tip, sample diluent, calibration QR code card and dilution vials</p>	<p>25 Tests</p>	<p>534020025</p>