



<p><b>FIAcheck™</b> <b>TT4</b></p>	<p><b>Fluorescence Immunoassay for quantitative determination of TT4 (Total Thyroxine) in human serum and plasma</b></p>
<p><b>Application</b></p>	<p>FIAcheck TT4 Fluorescence Immunoassay is intended for the in-vitro quantitative measurement of TT4 (Total Thyroxine) in human serum and plasma.</p>
<p><b>Principle</b></p>	<p>Rapid, lateral flow , quantitative immuno-fluorescence assay</p>
<p><b>Measuring range</b></p>	<p>10-350nmol/L</p>
<p><b>Reagents</b></p>	<p>Membrane assembly coated with</p> <ul style="list-style-type: none"> <li>• T4-KLH-Bovine IgG complex on test line</li> <li>• Sheep anti chicken IgY in control line</li> <li>• Mixture of TT4 mouse monoclonal antibody and chicken IgY labeled with fluorescent microspheres on the binding pad</li> </ul> <p>Sample diluent</p>
<p><b>Test duration</b></p>	<p>10 mins + 10 mins</p>
<p><b>Storage / Stability</b></p>	<p>4-30°C, 12 months</p>

Presentation	Pack	Cat.No.
<p><b>FIAcheck TT4</b> kit contains an individually pouched device and special tip, sample diluent, calibration QR code card and dilution vials</p>	<p>25 Tests</p>	<p>534020025</p>