

FIAcheckTM **1.0**



KEY FEATURES

- Europium labeled reagents for high sensitivity and specificity
- Faster turn around time
- Simple and convenient with one and two step reagent procedure
- Information management through direct LIS/HIS connection
- Calibration through secure QR code

TECHNICAL SPECIFICATIONS

Test Principle	Time Resolved Fluorescence Immunoassay
Measurement	Quantitative
Individual test speed	15s
Display	7-inch LCD screen
Connectivity	RS-232, RJ45 and Bluetooth
Printer	In-built thermal printer
Battery	In-built battery
Data transfer	USB port for data transfer and Software update
Dimension (LxWxH)	225mmx152mmx105mm
Weight	1.5Kg
Power supply	AC 100-240V, 50-60Hz

STANDARD EQUIPMENT LIST

Sr. No.	Components	Quantity
1	Power Adapter	1 pc
2	Qualified Certificate	1 pc
3	User Manual	1 pc
4	Quality control cards	2 pc
5	Typing paper	1 pc
6	Sample tube/Sample diluent holding rack	1 pc