



Dengucheck™ - NS 1

Rapid test for the detection of Dengue NS 1 antigen in human serum/plasma

Intended Use: _____

Dengucheck™ - NS 1 is a rapid, qualitative immunochromatographic test system for the detection of Dengue NS 1 (Dengue Non-Structural Protein-1) antigen in human serum or plasma. The test can be used as a screening test in early phase of Dengue viral infection in conjunction with other criteria. For professional use only.

Dengucheck™ - NS 1 kit components:

Device	Membrane assembly for Dengue NS-1 antigen is predisposed with: <ul style="list-style-type: none"> ● Agglutinating sera for dengue NS 1 antigen conjugated with colloidal gold ● Agglutinating sera for dengue NS 1 at the test window 'T' ● Agglutinating sera for mouse globulin at the control window 'C'
Dropper	Disposable plastic sample applicator
Other Accessories	Desiccant pouch and Package Insert

Performance Characteristics:

Sensitivity: 100%	Specificity: 100%
--------------------------	--------------------------

FEATURES	BENEFITS
Test system designed to detect dengue NS-1 antigen	Detects dengue infection on the first day of appearance of clinical symptoms.
No sample running buffer required	Easy, single step, user friendly procedure
Sensitivity 100% and specificity is 100%	Good sensitivity and specificity
Pack size 10 test with a shelf life of 24 months	Convenient for small laboratories also
Storage at 4-30°C	Well suited for on-site diagnosis, flexibility of remote distribution.

Storage / Stability	Temperature	Duration
Unopened kit/pouch	4-30°C	24 Months
Opened pouch	Room Temperature	Use immediately

Available Pack Sizes
10 Tests