



Falcivax™
Rapid test for the detection of Malaria Pv / Pf
(With alcohol swabs and blood lancets)

Intended Use: _____

Falcivax™ is a rapid test for the detection of *P. falciparum* specific histidine rich protein-2 (Pf HRP-2) and *P. vivax* specific pLDH in whole blood.

Falcivax™ kit components:

Device	Membrane assembly pre-dispensed with Agglutinating Sera for HRP-2 - colloidal gold conjugate, Agglutinating Sera for <i>P. vivax</i> specific pLDH - colloidal gold conjugate, rabbit globulin colloidal gold conjugate, Agglutinating Sera for HRP-2, Agglutinating Sera for Malaria specific pLDH and Agglutinating Sera for Rabbit globulin at the respective regions.
Buffer	Ready to use buffer in a dropper bottle
Sample applicator	Disposable inverted cup sample applicator
Other Accessories	Desiccant pouch, Package insert and Pictorial instructions for use

Performance Characteristics:

Sensitivity: 100% correlation with microscopy	Specificity: 100% to <i>P. falciparum</i> HRP-2 and <i>P. vivax</i> specific pLDH
--	--

FEATURES	BENEFITS
Detection and differentiation of <i>P. falciparum</i> and <i>P. vivax</i>	True speciation of the "Big Two" i.e. <i>P. falciparum</i> and <i>P. vivax</i>
Uses monoclonal antibodies against Pf. HRP-2 and Pv. pLDH	Sensitivity and specificity in excess of 98% for <i>P. falciparum</i> and <i>P. vivax</i>
Employs Heterophile Binding Reagent (HBR)	Tests results are interference protected
Uses whole blood from finger prick or primary tube, as a sample	Well suited for field conditions, on-site diagnosis, pathological labs and reference labs, alike.
Composite kit	No additional material required for testing
Easy to use	Simple 20 minutes assay, storage at 1-40° C

Storage / Stability	Temperature	Duration
Unopened kit/pouch	1°C - 40°C	24 Months
Opened pouch	Room Temperature	Use immediately

Available Pack Sizes	
10 Tests	25 Tests
50 Tests	100 Tests