



Insight CRP

Rapid Immunochromatographic assay for the detection of C-reactive protein in human serum.

Intended Use:

Insight CRP, a rapid, self-performing, Immunochromatographic assay for the detection of C-reactive protein in human serum.

Insight CRP Kit components:

Device	Membrane test assembly impregnated with colloidal gold conjugated to agglutinating sera for CRP and rabbit globulin, agglutinating sera for CRP and agglutinating sera for rabbit globulin at the respective regions.
Buffer	Sample running buffer
Sample Applicator	Disposable plastic sample applicator
Other Accessories	Desiccant pouch and Package Insert

Performance Characteristics:

Sensitivity: 0.6 mg/dl	Specificity: Anti-CRP antibodies specific to CRP
-------------------------------	---

FEATURES	BENEFITS
Cut-off at 6mg/L(0.6mg/dl)	Correct classification of positive and negative
Sample dilution (1:100)	Eliminate complement interference and risk of prozoning
Five minutes running time	Rapid turnaround time, quick diagnosis
Storage at 4-30°C	Flexibility of remote distribution, storage and on the spot testing even in physician's office and field settings

Storage / Stability	Temperature	Duration
Unopened kit	4 – 30°C	18 Months
Opened kit	4 – 30°C	Use immediately

Available Pack Sizes	
10 Tests	25 Tests