

**Insight DOA Panel 6.1****Rapid competitive immunochromatographic assay for the simultaneous detection of tetrahydrocannabinol, amphetamine, barbiturates, benzodiazepine, cocaine and opiates in human urine****Intended Use:** _____

Insight DOA panel 6.1, a rapid, self performing, qualitative, immunochromatographic assay for the simultaneous detection of tetrahydrocannabinol, amphetamine, barbiturates, benzodiazepine, cocaine and opiates in human urine.

Insight DOA 6.1 kit components:

Device	Membrane test assembly impregnated with colloidal gold conjugated to specific anti drug antibody (monoclonal) and mouse IgG, specific drug conjugated to BSA and goat antimouse IgG at the respective region
Dropper	Disposable plastic sample dropper
Other Accessories	Desiccant pouch and Package Insert

Sensitivity

THC: 50 ng/ml	BZO: 300 ng/ml	BAR: 300 ng/ml
AMP: 1000 ng/ml	COC: 300 ng/ml	OPI: 300 ng/ml

FEATURES	BENEFITS
High Efficiency	Directly detects the presence of the drug molecule in concentrations as low as 50 ng/ml or \geq 1000 ng/ml
Simple and Convenient	Ideal for Point of care testing
Urine used as a sample	Ensures that the sample can be collected easily without the need for a professional
Five minutes running time	Rapid turn around time, quick diagnosis
Storage at 4 – 30°C	Flexibility of remote distribution, storage and on the spot testing even in physician's office and field settings

Storage / Stability	Temperature	Duration
Unopened kit	4 – 30°C	18 Months
Opened kit	4 – 30°C	Use immediately

Available Pack Sizes
10 Tests