



Insight HBeAg

Rapid immunochromatographic assay for detection of HBeAg in human serum

Intended Use: _____

Insight HBeAg, a rapid, self performing, immunochromatographic assay that is used to detect the presence of HBeAg in human serum, thereby giving an assessment of the infectivity and the patient's status.

Insight HBeAg kit components:

Device	Membrane test assembly impregnated with colloidal gold conjugated to anti HBeAg monoclonal and mouse IgG antibodies, anti HBeAg (monoclonal) and goat anti mouse IgG at the respective regions
Dropper	Disposable plastic sample dropper
Other Accessories	Desiccant pouch and Package Insert

FEATURES	BENEFITS
Reliable Performance	Good correlation with commercial ELISA
Simple and Convenient	Rapid one step assay No need of expertise and automation
Fifteen minutes running time	Rapid turn around time, quick diagnosis
Storage at 4 – 30°C	Flexibility of remote distribution, storage and on the spot testing even in physician's office and field settings

Storage / Stability	Temperature	Duration
Unopened kit	4 – 30°C	18 Months
Opened pouch	4 – 30°C	Use immediately

Available Pack Sizes
10 Tests