



Insight PCP

Rapid competitive immunochromatographic assay for the detection of Phencyclidine in human urine

Intended Use: _____

Insight PCP, a rapid, self performing, qualitative, immunochromatographic assay which detects the presence of **Phencyclidine** in human urine specimen.

Insight PCP kit components:

Device	Membrane test assembly impregnated with colloidal gold conjugated to anti-phencyclidine antibody monoclonal and mouse IgG, phencyclidine conjugated to BSA and goat anti mouse IgG at the respective regions
Dropper	Disposable plastic sample dropper
Other Accessories	Desiccant pouch and Package Insert

Performance Characteristics:

Sensitivity: Detects the presence of phencyclidine in human urine specimens, qualitatively, at concentrations as low as 25 ng/ml.

FEATURES	BENEFITS
High Accuracy- well calibrated with toxicology controls	Directly detects the presence of the drug molecule in concentrations as low as 25 ng/ml
Simple and Convenient	Ideal for Point of care testing
Urine used as a sample	Ensures that the sample can be collected easily without the need for a professional settings
Five minutes running time	Rapid turn around time, quick diagnosis
Storage at 4 – 30°C	Flexibility of remote distribution, storage and on the spot testing even in physician's office and field settings

Storage / Stability	Temperature	Duration
Unopened kit	4 – 30°C	18 Months
Opened device	4 – 30°C	Use immediately

Available Pack Sizes
10 Tests