



## Tina™ hCG

### One step pregnancy test (Dipstick)

Intended Use: \_\_\_\_\_

Tina™ hCG is a rapid, qualitative two site sandwich immunoassay to detect presence of hCG in urine.

Tina™ hCG Kit components:

Dipstick	Membrane assembly pre-dispensed with anti-hCG antiserum-colloidal gold conjugate, rabbit IgG colloidal gold conjugate and anti-hCG antiserum and anti-rabbit antiserum at the respective regions
Other Accessories	Package Insert

Performance Characteristics:

Sensitivity: 100%	Specificity: 100% specific to -hCG
-------------------	------------------------------------

FEATURES	BENEFITS
Sensitivity calibrated to detect ~ 25 mIU/ml of hCG	Ideal for early pregnancy detection Exclude normal non pregnant individuals Facilitates safe and confident reporting of results
Simple and Accurate	One step test without requirement of skilled personnel Correlates very well with other commercially reliable kits
Highly specific	100% specificity for -hCG No cross reactivity with LH and FSH
Built in control	Automated validation of the test performance
Composite kit	30 individually pouched dipsticks in a container ensuring moisture-free protection
Results within 5 minutes	Quick and rapid
Multilingual test instructions	User friendly
Storage at 4-30°C	Facilitates remote distribution and on the spot testing

Storage/Stability	Temperature	Duration
Unopened kit	4-8°C	24 Months

Available Pack Sizes
30 Dipsticks