



<h1 style="color: green;">Cancheck-AFP</h1>	<h2>Rapid Test for AFP in human serum / plasma / whole blood</h2>
Application	For the detection of Alpha feto protein (AFP) levels in human serum/plasma/whole blood.
Principle	Self performing, semi-quantitative two-site sandwich immunoassay , lateral flow device
Sensitivity	100%
Specificity	99.45%
Reagents	Membrane assembly predispensed with <ul style="list-style-type: none"> • Agglutinating sera for AFP - colloidal gold conjugate and rabbit IgG colloidal gold conjugate • Agglutinating sera for AFP at 'Test' region (T) • Agglutinating sera for rabbit globulin at 'Control' Region (C)
Test duration	Five minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
Aluminium foil pouches containing CANCHECK-AFP membrane test device, disposable plastic dropper and sample running buffer	10 T	50900010

