

CANCHECK[®] - FOBT	Rapid Test for detection of faecal Ocult blood in human faeces
Application	For detection of faecal occult blood concentration in human faeces.
Principle	Self performing, semi-quantitative two-site sandwich immunoassay, lateral flow device
Sensitivity	Detection limit is upto 200ug/l of faeces suspension
Specificity	Highly specific to Human hemoglobin
Reagents	Membrane assembly predisposed with <ul style="list-style-type: none"> • Monoclonal Anti-human haemoglobin colloidal gold conjugate and rabbit IgG colloidal gold conjugate • Monoclonal Anti-human haemoglobin antibody and anti rabbit antiserum at the respective regions
Test duration	Five minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
Individual pouches each containing device and the sample dilution buffer	10 T	507040010

