



<b>Cancheck-PSA</b>	<b>Rapid Test for PSA in human serum/plasma/whole blood</b>
<b>Application</b>	For the detection of Prostate specific antigen (PSA) levels in human serum/plasma/whole blood.
<b>Principle</b>	Self performing, semi-quantitative two-site sandwich immunoassay , lateral flow device
<b>Sensitivity</b>	100%
<b>Specificity</b>	100%
<b>Reagents</b>	<p>Membrane assembly predispensed with</p> <ul style="list-style-type: none"> <li>• Monoclonal Anti-PSA antibody- colloidal gold conjugate and rabbit IgG colloidal gold conjugate</li> <li>• Specific Monoclonal Anti-PSA antibody at 'Test' region (T)</li> <li>• 10 ng/ml calibrated anti rabbit IgG antiserum at ` Reference region (R)</li> <li>• Anti-rabbit IgG antiserum at 'Control' region (C)</li> </ul>
<b>Test duration</b>	Five minutes assay
<b>Storage / Stability</b>	4-30°C, 2 years

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
Aluminium foil pouches containing <b>CANCHECK-PSA</b> membrane test device, disposable plastic dropper and sample running buffer	10 T	504010010