



	Rapid One Step Test for HCG
Application	For the detection of hCG in urine during pregnancy
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	~ 10 mIU/ml of hCG
Specificity	100% specific to hCG
Reagents	Membrane assembly predispensed with: <ul style="list-style-type: none"> - Agglutinating sera for hCG-colloidal gold conjugate, rabbit IgG colloidal gold conjugate - Agglutinating sera for hCG at 'Test' region - Agglutinating sera for rabbit globulin at 'Control' region
Test duration	Five minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
CLUE® kit contains individually pouched device with sample dropper	1 Test	301010001
	10 Tests	301010010
	25 Tests	301010025
	50 Tests	301010050
	100 Tests	301010100