



	<b>Rapid One Step Test for HCG</b>
<b>Application</b>	For the detection of hCG in urine during pregnancy
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow dipstick
<b>Sensitivity</b>	~ 10 mIU/ml of hCG
<b>Specificity</b>	100% specific to hCG
<b>Reagents</b>	Membrane assembly predisposed with: <ul style="list-style-type: none"> <li>– Agglutinating sera for hCG-colloidal gold conjugate, rabbit IgG colloidal gold conjugate</li> <li>– Agglutinating sera for hCG at 'Test' region</li> <li>– Agglutinating sera for rabbit globulin at 'Control' region</li> </ul>
<b>Test duration</b>	Five minutes assay
<b>Storage / Stability</b>	4-30°C, 2 years

<b>Presentation</b>	<b>Pack</b>	<b>Cat.No.</b>
<b>CLUE®</b> kit contains individually pouched dipstick	1 Test	301020001
	10 Tests	301020010
	25 Tests	301020025
	50 Tests	301020050
	100 Tests	301020100