



Coviscreen™	Rapid Double Antigen COVID-19 Screening test
Application	Detects Total antibodies (IgM+IgA+IgG) antibodies to SARS-CoV-2 virus in human serum/plasma/whole blood
Principle	A unique, self performing double antigen immunoassay on a membrane based on the principle of immunochromatography
Reagents	Membrane assembly predisposed with <ol style="list-style-type: none"> 1. SARS-CoV-2 virus specific recombinant antigen-indicator colloidal particles, streptavidin gold conjugate 2. SARS-CoV-2 virus specific recombinant antigen-capture at test region 'T' 3. Biotinylated BSA at the control region 'C'
Test duration	Twenty minutes assay
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Coviscreen™ kit contains individually pouched devices, sample droppers and sample running buffer	1 Test	502080001
	10 Tests	502080010
	25 Tests	502080025
	100 Tests	502080100