



<h1 style="color: green;">CoviRAT</h1>	Rapid test for the detection of SARS-CoV-2 antigen
Application	Detects SARS-CoV-2 antigen in human Nasopharyngeal Swab Specimens
Principle	A unique, self performing rapid, qualitative immunoassay on a membrane based on the principle of immunochromatography
Reagents	Membrane assembly predisposed with <ul style="list-style-type: none"> • Mouse monoclonal anti- SARS-CoV-2 antibodies--indicator colloidal particles • Mouse monoclonal anti-SARS-CoV-2 antibodies- capture at test region 'T' • Goat anti-mouse IgG at control region 'C'
Test duration	15-30 minutes assay
Storage / Stability	4-30°C, 15 months

Presentation	Pack	Cat.No.
CoviRAT kit contains individually pouched devices, Sterile nasopharyngeal swabs, Extraction Buffer tubes, Nozzles for extraction buffer tubes, Stands for extraction buffer tubes and Package insert	25 Tests	502100025

