



	Rapid test for the detection of Dengue NS 1 antigen in human serum/plasma
Application	Rapid test for the detection of NS-1 antigen to Dengue virus in human serum/plasma
Principle	A unique two-site, self performing immunoassay on a membrane based on the principle of immunochromatography
Sensitivity	NS-1=100%
Specificity	NS-1 = 100%
Reagents	Membrane assembly for Dengue NS-1 antigen is predisposed with: <ul style="list-style-type: none"> • Agglutinating sera for dengue NS 1 - colloidal gold conjugate • Agglutinating sera for dengue NS 1 at the test window 'T' • Agglutinating sera for mouse globulin predisposed at the control window 'C'
Test duration	Fifteen minutes assay
Storage / Stability	4-30°C, 24 months

Presentation	Pack	Cat.No.
Dengucheck™ NS1 kit contain device for detection of Dengue NS-1 antigen to dengue virus and sample running buffer along with sample dispensing dropper	10 Tests	502030010