



Dengucheck WB

Rapid Immunochromatographic Test for antibodies to Dengue virus

Application	New generation rapid immunochromatographic test for serodiagnosis of IgM and IgG antibodies to Dengue virus
Principle	A unique two-site, self performing immunoassay on a membrane based on the principle of immunochromatography
Sensitivity	100%
Specificity	100%
Reagents	<p>Dengucheck-WB kit contains:</p> <p>A. Individual pouches, each containing</p> <ol style="list-style-type: none"> 1. Test Device: Membrane assembly predisposed with Dengue virus specific antigen colloidal gold conjugate, streptavidin gold conjugate, anti human IgM at test region M, anti human IgG at test region G and biotinylated BSA at control region C. 2. Desiccant pouch 3. Disposable plastic sample dropper <p>Buffer bottles containing sample dilution buffer</p> <p>B. Package insert.</p>
Test duration	Twenty minute assay
Storage / Stability	4-30°C, 24 months

Presentation	Pack	Cat.No.
Dengucheck WB kit contains individually pouched devices and sample dilution buffer	10 Tests	50200010
	25 tests	50200025