



	<h2>Rapid immunochromatography test for HIV and Syphilis</h2>
Application	Rapid combo test for simultaneous detection of HIV 1/2 antibody and anti-treponemal antibodies in human serum, plasma or whole blood
Sensitivity	100% for HIV and Syphilis
Specificity	100% for HIV and Syphilis
Reagents	Membrane assembly predisposed with: <ul style="list-style-type: none"> • Recombinant HIV 1/2 antigen - colloidal gold conjugate • Recombinant T. pallidum antigen colloidal gold conjugate • Rabbit IgG-colloidal gold conjugate • Recombinant antigen for HIV 1 and HIV 2 at test line 'H' • Recombinant antigens specific to T. pallidum at the test line 'S' • Goat anti rabbit IgG coated at control region
Test duration	30 minutes
Storage / Stability	4°C to 30°C, 2 years

Presentation	Pack	Cat.No.
Duoscreen HIV/SYP kit contains individually pouched device and sample running buffer	25T	502170025