



FalciVax	Rapid Test for Malaria Pv/Pf
Application	For the detection of <i>P. falciparum</i> specific histidine rich protein-2 (Pf .HRP-2) and <i>P. vivax</i> specific pLDH in whole blood.
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	100% correlation with microscopy
Specificity	100% to <i>P.falciparum</i> HRP-2 and <i>P.vivax</i> specific pLDH
Reagents	<p>Membrane assembly pre-dispensed with</p> <ul style="list-style-type: none"> • Agglutinating Sera for HRP-2 - colloidal gold conjugate, • Agglutinating Sera for <i>P. vivax</i> specific pLDH - colloidal gold conjugate, rabbit globulin colloidal gold conjugate, • Agglutinating Sera for HRP-2, • Agglutinating Sera for Malaria specific pLDH • Agglutinating Sera for Rabbit globulin at the respective regions. <p>Clearing Buffer</p>
Test duration	Twenty minutes assay
Storage / Stability	1-40°C, 24 months

Presentation	Pack	Cat.No.
FalciVax kit contains <ul style="list-style-type: none"> • Individually pouched device with sample applicator pipette • Clearing buffer 	Mono Test	503150001
	5 Tests	503150005
	10 Tests	503150010
	25 Tests	503150025
	50 Tests	503150050

	100 Tests	503150100
--	-----------	-----------

