



| | |
|----------------------------|---|
| FalciVax | Rapid Test for Malaria Pv/Pf |
| Application | For the detection of <i>P. falciparum</i> specific histidine rich protein-2 (Pf .HRP-2) and <i>P. vivax</i> specific pLDH in whole blood. |
| Principle | Self performing, two site sandwich immunoassay, lateral flow device |
| Sensitivity | 100% correlation with microscopy |
| Specificity | 100% to <i>P.falciparum</i> HRP-2 and <i>P.vivax</i> specific pLDH |
| Reagents | <p>Membrane assembly pre-dispensed with</p> <ul style="list-style-type: none"> • Agglutinating Sera for HRP-2 - colloidal gold conjugate, • Agglutinating Sera for <i>P. vivax</i> specific pLDH - colloidal gold conjugate, rabbit globulin colloidal gold conjugate, • Agglutinating Sera for HRP-2, • Agglutinating Sera for Malaria specific pLDH and • Agglutinating Sera for Rabbit globulin at the respective regions. <p>Clearing Buffer</p> |
| Test duration | Twenty minutes assay |
| Storage / Stability | 1-40°C, 24 months |

| Presentation | Pack | Cat.No. |
|---|-------------|----------------|
| FalciVax kit contains <ul style="list-style-type: none"> • Individually pouched device with sample applicator pipette • Clearing buffer • Alcohol Swabs and blood lancets | Mono Test | 503010001 |
| | 10 Tests | 503010010 |
| | 25 Tests | 503010025 |
| | 50 Tests | 503010050 |
| | 100 Tests | 503010100 |

