



<h1>Ferticheck-LH</h1>	<h2>Rapid One Step Test for human Luteinising hormone in urine</h2>
<b>Application</b>	<p>For the detection of human Luteinising hormone (hLH) concentration in urine.</p> <p>A 7-day ovulation test designed to predict the hLH surge associated with the time of ovulation and best probable time to plan conception.</p>
<b>Principle</b>	<p>Self performing, semi-quantitative two-site sandwich immunoassay , lateral flow device</p>
<b>Sensitivity</b>	<p>≥ 40 mIU/ml of hLH in urine</p>
<b>Specificity</b>	<p>No cross reactivity with human follicle stimulating hormone (hFSH) and human thyroid stimulating hormone (hTSH)</p>
<b>Reagents</b>	<ul style="list-style-type: none"> <li>• Membrane assembly predisposed with             <ul style="list-style-type: none"> <li>- Monoclonal Anti-hLH antibody- colloidal gold conjugate and rabbit IgG colloidal gold conjugate</li> <li>- anti alpha hCG antibody at 'Test' region (T)</li> <li>- Anti-rabbit antiserum at 'Control' Region (C)</li> </ul> </li> </ul>
<b>Test duration</b>	<p>Five minutes assay</p>
<b>Storage / Stability</b>	<p>4-30°C, 2 years</p>

Presentation	Pack	Cat.No.
<p>Aluminium foil pouches containing <b>FERTICHECK-LH</b> membrane test assembly and disposable plastic dropper</p>	<p>5 Tests</p>	<p>505010005</p>