



<h1 style="color: green;">Ferticheck-LH</h1>		<h2>Rapid One Step Test for human Luteinising hormone in urine</h2>
Application	For the detection of human Luteinising hormone (hLH) concentration in urine. A 7-day ovulation test designed to predict the hLH surge associated with the time of ovulation and best probable time to plan conception.	
Principle	Self performing, semi-quantitative two-site sandwich immunoassay , lateral flow device	
Sensitivity	≥ 40 mIU/ml of hLH in urine	
Specificity	No cross reactivity with human follicle stimulating hormone (hFSH) and human thyroid stimulating hormone (hTSH)	
Reagents	<ul style="list-style-type: none"> • Membrane assembly predispensed with <ul style="list-style-type: none"> - Monoclonal Anti-hLH antibody-colloidal gold conjugate and rabbit IgG colloidal gold conjugate - anti alpha hCG antibody at 'Test' region (T) - Anti-rabbit antiserum at 'Control' Region (C) 	
Test duration	Five minutes assay	
Storage / Stability	4-30°C, 2 years	

Presentation	Pack	Cat.No.
Aluminium foil pouches containing FERTICHECK-LH membrane test assembly and disposable plastic dropper	5 Tests	51100005