



	<b>Rapid test for the detection of total antibodies to HCV in serum/plasma/whole blood</b>
<b>Application</b>	For the detection of total antibodies to HCV in serum/plasma/whole blood
<b>Principle</b>	Advanced, fourth generation, rapid, self performing, qualitative, two site, sandwich lateral flow immunochromatographic assay
<b>Sensitivity</b>	100%
<b>Specificity</b>	100%
<b>Reagents</b>	<ul style="list-style-type: none"> <li>• Membrane assembly predisposed with multi-epitope HCV Recombinant Peptide Antigen-colloidal gold conjugate, rabbit IgG colloidal gold conjugate</li> <li>• Multi-epitope HCV specific Recombinant Peptide Antigen at region 'T'</li> <li>• Agglutinating sera for rabbit globulin at region 'C'</li> </ul>
<b>Test duration</b>	Fifteen minutes
<b>Storage / Stability</b>	4-30°C, 2 years

Presentation	Pack	Cat.No.
<b>FLAVICHECK-HCV WB</b> kit contains individually pouched devices and sample running buffer	12 x MONO TEST	40701012
	1 TEST	40701001
	10 TESTS	40701010
	25 TESTS	40701025
	50 TESTS	40701050

