



	Rapid immunochromatography test for HCV
Application	For the detection of total antibodies (IgM, IgG, IgA) to HCV in Human serum / plasma
Principle	Self-performing, two-site sandwich immunoassay, lateral flow device
Sensitivity	100%
Specificity	99.61%
Reagents	Membrane assembly predisposed with: <ul style="list-style-type: none"> • Multi-epitope HCV recombinant antigen conjugated to colloidal gold, rabbit IgG conjugated to colloidal gold. • Multi-epitope recombinant antigen at 'Test' region. Agglutinating sera for rabbit globulin at 'Control' region.
Test duration	15 minutes assay
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
FLAVISCREEN-HCV PLUS test device and sample running buffer	10 TESTS	402170010
	50 TESTS	402170050