



	Rapid One Step Test for HCG
Application	For the detection of hCG in urine/serum during pregnancy
Principle	Self performing, double antigen sandwich, immunoassay, lateral flow device
Sensitivity	~ 10 mIU/ml of hCG
Specificity	100% specific to b hCG
Reagents	<ul style="list-style-type: none"> • Membrane assembly predisposed with <ul style="list-style-type: none"> - Monoclonal anti-hCG antisera-colloidal gold conjugate, rabbit IgG colloidal gold conjugate - Monoclonal anti-hCG antiserum at the 'Test region' - Anti-rabbit anti serum at the 'Control region'
Test duration	Urine 5 mins. Serum 15 mins.
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
GRAVICHECK® kit contains individually pouched device with dropper	25 Tests	40201025

