



	Rapid One Step Test for HCG
Application	For the detection of hCG in urine/serum during pregnancy
Principle	Self performing, two site sandwich, qualitative, immunoassay, lateral flow dipstick
Sensitivity	~ 10 mIU/ml of hCG
Specificity	100% specific to hCG
Reagents	<ul style="list-style-type: none"> • Membrane assembly predisposed with <ul style="list-style-type: none"> - Monoclonal Anti-hCG antisera-colloidal gold conjugate, rabbit IgG colloidal gold conjugate - Monoclonal Anti-hCG antiserum at 'Test' region - Anti-rabbit antiserum at 'Control' region
Test duration	Urine 5 mins. Serum 15 mins.
Storage / Stability	4-30°C, 2 years

Presentation	Pack	Cat.No.
GRAVICHECK® kit contains individually pouched dipstick	25 Tests	40202025

