



<h1 style="color: green; margin: 0;">HelicoCheck™</h1>	<h2 style="margin: 0;">Rapid Test for detection of antibodies to <i>Helicobacter pylori</i> in human serum/plasma</h2>
<b>Application</b>	For the detection of antibodies to <i>Helicobacter pylori</i> in human serum or plasma.
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow device.
<b>Sensitivity</b>	92.3%
<b>Specificity</b>	90.5%
<b>Reagents</b>	<ul style="list-style-type: none"> <li>Membrane test assembly impregnated with colloidal gold conjugated to recombinant <i>Helicobacter pylori</i> antigen (tracer) and rabbit IgG</li> <li>Recombinant <i>Helicobacter pylori</i> antigen (capture) and anti rabbit antisera at the respective regions</li> </ul>
<b>Test duration</b>	Fifteen minutes assay
<b>Storage / Stability</b>	4-30°C, 2 years

Presentation	Pack	Cat.No.
<b>Helicocheck kit</b> contains individually pouched device with sample dispensing droppers	10T	501070010
	25T	518010025

