



	<h2>One Step Test for HBsAg</h2>
<b>Application</b>	For the detection of Hepatitis B surface antigen, A, a marker for Hepatitis B infections in serum/ plasma specimen.
<b>Principle</b>	Rapid, qualitative, two site sandwich immunoassay, lateral flow device
<b>Detection limit</b>	~ 0.5 ng/ml of HBsAg
<b>Sensitivity</b>	100%
<b>Specificity</b>	100%
<b>Reagents</b>	Membrane assembly predispensed with: <ul style="list-style-type: none"> <li>• Agglutinating Sera for HBsAg             <ul style="list-style-type: none"> <li>- colloidal gold conjugate, rabbit IgG</li> <li>colloidal gold conjugate</li> </ul> </li> <li>• Agglutinating Sera for HBsAg at 'Test' region</li> <li>• Agglutinating Sera for Rabbit globulin at 'Control' region</li> </ul>
<b>Test duration</b>	15 mins.
<b>Storage / Stability</b>	4-30°C, 2 years

Presentation	Pack	Cat.No.
<b>HEPAVIEW</b> kit contains individually pouched device with sample dropper	25 TESTS	402070025

