



	Rapid Immunochromatographic assay for the detection of Anti-Streptolysin O in human serum
Application	For the detection of Anti-Streptolysin O in human serum
Sensitivity	240 IU/ml
Reagents	<p>A. Each INSIGHT ASO kit contains individual pouches each containing:</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to SLO antigen and the streptavidin, agglutinating sera for human IgG and biotinylated BSA at the respective regions. 2. Desiccant pouch. 3. Sample applicator. 4. sample running buffer <p>B. Package insert.</p>
Test duration	Five minutes assay
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ ASO kit contains individually pouched device	10 Tests	108190010
	25 Tests	108190025

