



	<b>Rapid Competitive Immunochemical Test for detection of Benzodiazepine in urine</b>
<b>Application</b>	For the detection of presence of Benzodiazepine in human urine specimens qualitatively with concentrations as low as 300 ng/ml
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow device
<b>Reagents</b>	A. Each INSIGHT BZO kit contains individual pouches each containing a <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti benzodiazepine antibody monoclonal and mouse IgG, benzodiazepine conjugated to BSA and goat anti mouse IgG at the respective regions.</li> <li>2. Sample dropper.</li> <li>3. Desiccant pouch.</li> </ol> B. Package insert.

Presentation	Pack	Cat.No.
<b>Insight™ BZO kit contains individually pouched device</b>	10 Tests	10808010

