



	Rapid Competitive Immunochromatographic Test for detection of cocaine in urine
Application	For the detection of presence of cocaine in human urine specimens qualitatively with concentrations as low as 300 ng/ml
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Reagents	<p>A. Each INSIGHT COC kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti cocaine antibody monoclonal and mouse IgG, cocaine conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p>

Presentation	Pack	Cat.No.
Insight™ COC kit contains individually pouched device	10 Tests	10802010