



	Rapid Immunochromatographic assay for the detection of C-reactive protein in human serum
Application	For the detection of C-reactive protein in human serum
Sensitivity	6 mg/L (0.6 mg/dl)
Reagents	<p>A. Each INSIGHT CRP kit contains individual pouches each containing:</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to agglutinating sera for CRP and rabbit globulin, agglutinating sera for CRP and agglutinating sera for rabbit globulin at the respective regions. 2. Desiccant pouch. 3. Sample applicator. 4. sample running buffer <p>B. Package insert.</p>
Test duration	Five minutes assay
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ CRP kit contains individually pouched device	10 Tests	108170010
	25 Tests	108170025

