



<p>Insight DOA Panel 6.1 DEVICE</p>	<p>Rapid Competitive Immunochromatographic assay for the simultaneous detection of tetrahydrocannabinol, amphetamine, barbiturates, benzodiazepine, cocaine and opiates in human urine</p>
<p>Application</p>	<p>For the simultaneous detection of tetrahydrocannabinol, amphetamine, barbiturates, benzodiazepine, cocaine and opiates in human urine</p>
<p>Principle</p>	<p>Self performing, two site sandwich immunoassay, lateral flow device</p>
<p>Reagents</p>	<p>A. Each INSIGHT MTD kit contains individual pouches containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to specific anti drug antibody (monoclonal) and mouse IgG, specific drug conjugated to BSA and goat antimouse IgG at the respective region. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p>
<p>Sensitivity</p>	<p>≥ 50 ng/ml of THC ≥ 1000 ng/ml of Ampethamine ≥ 300 ng/ml of BAR ≥ 300 ng/ml of BZO</p>

	≥ 300 ng/ml of COC ≥ 300 ng/ml of OPI
--	--

Presentation	Pack	Cat.No.
Insight™ DOA 6.1 kit contains individually pouched devices	10 Tests	10810010

