



	Rapid Competitive ImmunoChromatographic Test for detection of HBcAb
Application	For the detection of HBcAb in human serum
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	~ 4 ncu/ml
Reagents	<p>A. Each INSIGHT HBcAb kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti HBcAb monoclonal and mouse IgG antibodies, HBcAg, anti HBcAg (monoclonal) and goat anti mouse IgG at the respective regions. 2. Desiccant pouch. 3. Sample dropper. <p>B. Package insert.</p>
Test duration	15 mins.
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ HBcAb kit contains individually pouched device	10 Tests	10503010