



	Rapid Competitive Immunochromatographic Test for detection of HBeAb
Application	For the detection of HBeAb in human serum
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	~ 4.0 ncu/ml
Reagents	A. Each INSIGHT HBeAb kit contains individual pouches each containing a <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to HBeAg monoclonal antibody and mouse IgG, HBeAg, anti HBeAg (monoclonal) antibody and goat anti mouse IgG at the respective regions. 2. Desiccant pouch. 3. Sample dropper. B. Package insert.
Test duration	15 mins.
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ HBeAb kit contains individually pouched device	10 Tests	10503010

