



	<b>Rapid Immunochromatographic Test for detection of HBeAg</b>
<b>Application</b>	For the detection of HBeAg in human serum
<b>Principle</b>	Self performing, two site sandwich immunoassay, lateral flow device
<b>Sensitivity</b>	$\sim 1$ ncu/ml
<b>Reagents</b>	<p>A. Each INSIGHT HBeAg kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> <li>1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti HBeAg monoclonal and mouse IgG antibodies, anti HBeAg (monoclonal) and goat anti mouse IgG at the respective regions.</li> <li>2. Desiccant pouch.</li> <li>3. Sample dropper.</li> </ol> <p>B. Package insert.</p>
<b>Test duration</b>	15 mins.
<b>Storage / Stability</b>	4-30°C, 18 months

Presentation	Pack	Cat.No.
<b>Insight™ HBeAg kit contains individually pouched device</b>	10 Tests	10506010

