



	Rapid Immunochromatographic Test for detection of HBsAb
Application	For the detection of HBsAb in human serum
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Sensitivity	~ 30 mIU
Reagents	A. Each INSIGHT HBsAb kit contains individual pouches each containing a <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to HBsAg and mouse IgG antibodies, HBsAg antigen and goat anti mouse IgG at the respective regions. 2. Desiccant pouch. 3. Sample dropper. B. Package insert.
Test duration	15 mins.
Storage / Stability	4-30°C, 18 months

Presentation	Pack	Cat.No.
Insight™ HBsAb kit contains individually pouched device	10 Tests	10505010

