



	Rapid Competitive Immunochromatographic Test for detection of Ketamine in urine
Application	Detects the presence of Ketamine in human urine at cutoff concentration 250 ng/ml and above
Principle	Self performing, two site sandwich immunoassay, lateral flow device
Reagents	<p>A. Each INSIGHT KET kit contains individual pouches each containing a</p> <ol style="list-style-type: none"> 1. Device: Membrane test assembly impregnated with colloidal gold conjugated to anti ketamine antibody monoclonal and mouse IgG, ketamine conjugated to BSA and goat anti mouse IgG at the respective regions. 2. Sample dropper. 3. Desiccant pouch. <p>B. Package insert.</p>

Presentation	Pack	Cat.No.
Insight™ KET kit contains individually pouched device	10 Tests	10805010

